

thermoscientific



Hypersil GOLD HPLC Columns

Outstanding peak shape for your separations

 fisher scientific

ThermoFisher
SCIENTIFIC

Hypersil GOLD Columns

Designed for improved chromatography, Thermo Scientific™ Hypersil GOLD™ columns are the culmination of 40 years of experience in the product development and manufacturing of HPLC media and columns. The range and capabilities of this state-of-the-art family of columns, with numerous chemistries and a range of particle sizes and hardware formats meet the challenges of modern chromatography.

The highly pure Hypersil GOLD silica is manufactured, bonded and packed in ISO 9001:2008 accredited facilities, operating under strict protocols using robust procedures and extensive quality control testing. The manufacturing and bonding process creates an even surface with fewer silanols leading to reduced secondary interactions. This ensures consistent performance, column after column.

Hypersil GOLD HPLC columns are available in 12 different chemistries to optimize separations and maximize productivity. The extensive range of Hypersil GOLD columns offers chromatographers outstanding peak shape for reversed phase, ion exchange, HILIC or normal phase chromatography. With all 12 phases being available with 1.9 µm particle size, Hypersil GOLD columns offers chromatographers flexibility in choosing the correct column, whether they are using conventional or ultra-high pressure LC systems.




Improved Selectivity, Resolution and Productivity

Hypersil GOLD Outstanding Peak Shape Using Generic Gradients with C18 Selectivity	2–7
Hypersil GOLD C8 Enhanced Resolution, Efficiency, Sensitivity and Speed	8–9
Hypersil GOLD C4 Low Hydrophobicity Columns for Less Retention	10
Hypersil GOLD aQ Enhanced Retention and Resolution of Polar Analytes	11–12
Hypersil GOLD PFP Unique Selectivity with Perfluorinated Columns	13–14
Hypersil GOLD CN Cyano Columns for Reversed and Normal Phase Separations	15–16
Hypersil GOLD Phenyl Excellent Retention and Unique Selectivity for Aromatic Analytes	17–18
Hypersil GOLD Amino Highly Versatile Aminopropyl Stationary Phase	19
Hypersil GOLD AX Separation of Anionic Species and Polar Molecules	20
Hypersil GOLD SAX Quarternary Amine Strong Anion Exchange Column, Designed for Aqueous Mobile Phase	21
Hypersil GOLD Silica Excellent Peak Shape in Normal Phase Chromatography	22
Hypersil GOLD HILIC Enhanced Retention of Polar and Hydrophilic Analytes	23
Hypersil GOLD 1.9 μm Small Particles to Improve Speed and Efficiency	24–26
Hardware Solutions Hardware Solutions for High Throughput Screening, Capillary and Preparative Chromatography	27
Column Protection Extend Column Lifetime and Improve Performance	28
Ordering Information	29–32


Hypersil GOLD

Outstanding peak shape using generic gradients with C18 selectivity



Hypersil GOLD columns are based on highly pure silica and a novel proprietary derivatization and endcapping procedure using alkyl chain chemistry. This gives:


- Significant reduction in peak tailing while retaining C18 (USP L1) selectivity
- Excellent resolution, efficiency and sensitivity
- Confidence in the accuracy and quality of analytical data



Hypersil GOLD columns offer improved peak shape, even for basic analytes.

Enhanced Resolution


Robust assay development requires a clear definition of resolution expectations. Narrow symmetrical chromatographic peaks ensure that optimum resolution is achieved. Obtaining narrow peak widths is especially challenging for basic pharmaceutical compounds. The reduced silanol activity on Hypersil GOLD columns reduces tailing for basic analytes, thus improving resolution.



Hypersil GOLD columns provide excellent resolution between critical pairs, aiding separation of closely related species.


Improved Sensitivity

Outstanding peak shape results in greater sensitivity. When peaks exhibit tailing, peak height is reduced, therefore compromising the sensitivity of the analysis. The highly symmetrical peaks provided by Hypersil GOLD columns enhance peak height and allow for optimised peak integration calculations. This can be particularly critical when low concentrations of an analyte are present, for example in an impurity assay.



Reproducibility

Our Hypersil GOLD columns are exceptionally reproducible for reliable chromatography, column after column. This allows the user to be confident that assays developed with Hypersil GOLD columns will be robust and stable for the life of the assay, making them an ideal choice for new method development.




Column: Hypersil GOLD , 5 μ m, 150 \times 4.6 mm
 Mobile Phase: A: 0.1% ammonia pH 10.6
 B: Methanol + 0.1% ammonia
 Gradient: 5–100% B in 15 min
 Flow Rate: 1.0 mL/min
 Injection Volume: 10 μ L
 Detection: UV @ 254 nm
 Temperature: 30 $^{\circ}$ C

High pH stability assay (pH 10.6) of Hypersil GOLD columns.

pH Stability


Our Hypersil GOLD columns are well suited to extended pH applications and have been shown to produce robust assays at high pH. At low pH, excellent column stability and reproducibility are illustrated over 1500 injections at pH 1.8.



Stability of Hypersil GOLD columns at low pH. No loss of retention after 28 L of mobile phase in 19.5 days of analysis.

Pharmaceutical


Cepha antibiotics



Column: Hypersil GOLD, 5 μ m, 150 \times 4.6 mm
 Mobile Phase: A: 0.1% acetic acid
 B: Acetonitrile
 Gradient: 20–70% B in 10 mins
 Flow Rate: 1 mL/min
 Detection: UV @ 254 nm
 Temperature: 25 $^{\circ}$ C

1. Cefadroxil
2. Cefaclor
3. Cephalexin
4. Cephradine
5. Cefazolin

Cough/cold formulation




Column: Hypersil GOLD, 5 μ m, 150 \times 4.6 mm
 Mobile Phase: A: 20 mM ammonium formate at pH 3.0
 B: Methanol
 Gradient:

Time (min)	% B
0	10
5	10
10	70

 Flow Rate: 1.5 mL/min
 Detection: UV @ 270 nm
 Temperature: 25 $^{\circ}$ C

1. 4-Amino phenol
2. (chlorpheniramine) maleate
3. Phenylephrine
4. Acetaminophen
5. Saccharin
6. Impurity from 4-Amino phenol
7. 4-Nitro phenol
8. Chlorpheniramine

Anaesthetics




Column: Hypersil GOLD, 5 μ m, 150 \times 4.6 mm
 Mobile Phase: A: 0.05 M monopotassium phosphate pH 3
 B: Acetonitrile
 Isocratic: 50:50
 Flow Rate: 1.25 mL/min
 Detection: UV @ 220 nm
 Temperature: 25 $^{\circ}$ C

1. Lidocaine
2. Tetracaine
3. Benzocaine

Environmental


Polycyclic aromatic hydrocarbons



Column: Hypersil GOLD, 5 μ m, 150 \times 4.6 mm
 Mobile Phase: A: Methanol
 B: Water
 Isocratic: 75:25
 Flow Rate: 1 mL/min
 Detection: UV @ 269 nm
 Temperature: 25 $^{\circ}$ C

1. Naphthalene
 2. Fluorene
 3. Phenanthrene
 4. Anthracene
 5. Pyrene
 6. Chrysene


Banned aromatic amines



Column: Hypersil GOLD, 3 μ m, 150 \times 2.1 mm
 Mobile Phase: A: 25 mM ammonium acetate at pH 5
 B: Acetonitrile
 Gradient: 20–100% B in 10 min
 Flow Rate: 0.2 mL/min
 Detection: UV @ 254 nm
 Temperature: 40 $^{\circ}$ C

1. 2,4-Diaminotoluene
 2. 4,4'-Oxydianiline
 3. o-Toluidine
 4. 2-Methoxy-5-methylaniline
 5. 2,4,5-Trimethylaniline
 6. 4,4'-Methylene-bis(2-chloroaniline)
 7. Unknown

Endocrine disruptors




Column: Hypersil GOLD, 5 μ m, 150 \times 4.6 mm
 Mobile Phase: A: Water
 B: Acetonitrile
 Gradient: 25–70% B in 20 min
 Flow Rate: 1.5 mL/min
 Detection: UV @ 220 nm
 Temperature: 25 $^{\circ}$ C

1. Desethyl atrazine
 2. Estriol
 3. Simazine
 4. Atrazine
 5. Diuron
 6. Bisphenol A
 7. Estrone

Toxicology


Testosterones



Column: Hypersil GOLD, 5 μ m, 150 \times 4.6 mm
 Mobile Phase: A: Water
 B: Acetonitrile
 Isocratic: 43:57
 Flow Rate: 1 mL/min
 Detection: UV @ 254 nm
 Temperature: 25 $^{\circ}$ C

- 1. 11-Ketotestosterone
- 2. 19-Nortestosterone (nandrolone)
- 3. Testosterone
- 4. Epitestosterone


Chlorpromazine



Column: Hypersil GOLD, 5 μ m, 50 \times 2.1 mm
 Mobile Phase: A: 0.1% formic acid
 B: Acetonitrile + 0.1% formic acid
 Gradient: 15–80% B in 5 min
 Flow Rate: 1 mL/min
 Detection: UV @ 254 nm
 Temperature: 30 $^{\circ}$ C

- 1. Chlorpromazine

Tricyclic antidepressants




Column: Hypersil GOLD, 5 μ m, 150 \times 4.6 mm
 Mobile Phase: A: 0.1% formic acid
 B: Acetonitrile + 0.1% formic acid
 Gradient: 30–50% B in 15 min
 Flow Rate: 1 mL/min
 Detection: UV @ 254 nm
 Temperature: 30 $^{\circ}$ C
 Concentration: 2.5 ng/ μ L

- 1. Doxepin
- 2. Protriptyline
- 3. Imipramine
- 4. Nortriptyline
- 5. Amitriptyline
- 6. Trimipramine

Food Safety


Energy drink additives



Column: Hypersil GOLD, 5 μ m, 150 \times 4.6 mm
 Mobile Phase: A: 10 mM ammonium acetate at pH 5.0
 B: Methanol
 Gradient: 30–45% B in 10 min
 Flow Rate: 1 mL/min
 Detection: UV @ 230 nm
 Temperature: 25 $^{\circ}$ C

1. Acesulfame
2. Saccharin
3. Caffeine
4. Benzoic acid
5. Sorbic acid
6. Aspartame


Coumaric acids



Column: Hypersil GOLD, 5 μ m, 150 \times 4.6 mm
 Mobile Phase: A: 0.1% formic acid
 B: Acetonitrile
 Isocratic: 70:30
 Flow Rate: 1 mL/min
 Detection: UV @ 270 nm
 Temperature: 40 $^{\circ}$ C

1. Uracil
2. p-Coumaric Acid
3. m-Coumaric Acid
4. o-Coumaric Acid

Tocopherols




Column: Hypersil GOLD, 5 μ m, 150 \times 4.6 mm
 Mobile Phase: A: Water
 B: Methanol
 Isocratic: 5:95
 Flow Rate: 1 mL/min
 Detection: -ESI
 Temperature: 30 $^{\circ}$ C

1. δ -Tocopherol
2. γ -Tocopherol
3. α -Tocopherol

Hypersil GOLD C8

Enhanced resolution, efficiency, sensitivity and speed



Particle Size 1.9 μm , 3 μm , 5 μm **USP** L7

- Analytes of medium hydrophobicity
- When a less hydrophobic phase is required to obtain adequate retention

Similar Selectivity but Less Retention than C18

Hypersil GOLD C8 media provides similar selectivity to C18 with a predictable elution order, but less retention. This feature is particularly useful where lower hydrophobicity is needed in order to successfully retain compounds of interest. Hypersil GOLD C8 columns are recommended for analytes of medium hydrophobicity or when a less hydrophobic phase is required to obtain adequate retention.

Faster Separations


Hypersil GOLD C8 columns can provide improved throughput of analysis over that of a C18 alkyl chain chemistry. Hydrophobic interactions are reduced, allowing compounds to elute quicker from the column.

Excellent Peak Shapes with High Efficiency and Outstanding Sensitivity

Hypersil GOLD C8 columns provide very symmetrical peak shapes while also improving capabilities such as speed of analysis, efficiency and sensitivity.




Food Safety

 β -Carotene

Column: Hypersil GOLD C8, 5 μ m, 150 \times 4.6 mm
 Mobile Phase: Methanol
 Flow Rate: 1.5 mL/min
 Detection: UV @ 450 nm
 Temperature: 25 $^{\circ}$ C

1. Lutein
 2. Lycopene
 3. β -Carotene

Fatty acids




Column: Hypersil GOLD C8, 5 μ m, 150 \times 4.6 mm
 Mobile Phase: A: 0.1% formic acid
 B: Acetonitrile
 Isocratic: 15:85
 Flow Rate: 1 mL/min
 Detection: UV @ 200 nm
 Temperature: 25 $^{\circ}$ C

1. Linolenic Acid
 2. Linoleic Acid

Environmental

Triazines and uron herbicides




Column: Hypersil GOLD C8, 5 μ m, 150 \times 4.6 mm
 Mobile Phase: A: Water
 B: Acetonitrile
 Gradient:

Time (min)	% B
0	20
15	23
25	75

Flow Rate: 1.5 mL/min
 Detection: UV @ 240 nm
 Temperature: 25 $^{\circ}$ C

1. Simazine
 2. Monuron
 3. Chlorotoluron
 4. Atrazine
 5. Diuron
 6. Propazine
 7. Linuron

Phthalates




Column: Hypersil GOLD C8, 5 μ m, 150 \times 4.6 mm
 Mobile Phase: A: Water
 B: Acetonitrile
 Gradient: 60–90% B in 10 min; hold 10 min
 Flow Rate: 1 mL/min
 Detection: UV @ 254 nm
 Temperature: 25 $^{\circ}$ C

1. Dimethyl phthalate
 2. Diethyl phthalate
 3. Dipropyl phthalate
 4. Diisopropyl phthalate
 5. Di-n-butyl phthalate
 6. Di-n-octyl phthalate

Hypersil GOLD C4

Low hydrophobicity columns for less retention



- Analytes with high hydrophobicity
- When a less hydrophobic phase is required to obtain adequate retention

Lower Hydrophobicity for Faster Separations


Hypersil GOLD C4 columns provide similar selectivity to C18 and C8 columns but with less retention. The shorter chain length and lower hydrophobic character make C4 a particularly useful stationary phase for the retention and separation of hydrophobic polypeptides and proteins.

Excellent Peak Shape, Showing High Efficiency and Outstanding Sensitivity

Based on the same highly pure silica, Hypersil GOLD C4 columns deliver excellent peak shape. For high speed, high efficiency separations, Hypersil GOLD C4 columns are available with 1.9 µm particle size.

Pharmaceutical

Parabens




Column: Hypersil GOLD C4, 5 µm, 150 × 4.6 mm
Mobile Phase: Water/acetonitrile (50:50)
Flow Rate: 1.0 mL/min
Temperature: 25 °C
Detection: 214 nm
Injection volume: 10 µL

1. Methylparaben
2. Ethylparaben
3. Propylparaben
4. Butylparaben

Food safety

Fatty Acids




Column: Hypersil GOLD C4, 1.9 µm, 100 × 2.1 mm
Mobile Phase: Water/acetonitrile (20:80)
Flow Rate: 0.55 mL/min
Temperature: 30 °C
Detection: 200 nm
Injection volume: 1 µL

1. Linolenic acid
2. Linoleic acid

Hypersil GOLD aQ

Enhanced retention and resolution of polar analytes



- Analysis of water soluble vitamins and organic acids
- Use with highly aqueous mobile phase

Retention and Resolution of Polar Analytes

Because Hypersil GOLD aQ columns are packed with a polar endcapped C18 phase, they offer superior retention of polar compounds. Dispersive interactions are the primary mechanism of retention with alkyl chain bonded phases. The polar functional group used to endcap Hypersil GOLD aQ media provides an additional controlled interaction mechanism by which polar compounds can be retained and resolved. The resulting optimized peak shape provides excellent resolution sensitivity and efficiency, making Hypersil GOLD aQ columns ideal for the quantitative analysis of trace levels of polar analytes.

Polar Endcapped C18 Stationary Phase for Alternative Selectivity

The additional interaction mechanism often provides selectivity differences over the traditional alkyl chain chemistries, and offers a solution for the separation of polar compounds which exhibit insufficient retention on pure alkyl chain phases under typical reversed phase mobile phase conditions.

Ideal for Highly Aqueous Mobile Phases


The wettability of reversed phase media can be increased by the introduction of polar functional groups. The polar endcapping of Hypersil GOLD aQ media also makes it usable in 100% aqueous mobile phases without the risk of loss of performance or poor stability.

Excellent Peak Shapes

Hypersil GOLD aQ silica ensures optimized peak shape, resolution, sensitivity and efficiency. Hypersil GOLD aQ columns provide only controlled secondary interactions to ensure excellent peak shape for all analyte types, making them ideal for the quantitative analysis of trace levels of polar analytes.

Food Safety

Water soluble vitamins




Column: Hypersil GOLD aQ, 5 μ m, 150 \times 4.6 mm
 Mobile Phase: A: 50 mM monopotassium phosphate pH 3.5
 B: Methanol
 Gradient: 0–100% B in 15 min
 Flow Rate: 1 mL/min
 Detection: UV @ 205 nm

1. Vitamin B1 (thiamine)
2. Vitamin B6 (pyridoxine)
3. Vitamin B3 (nicotinamide)
4. Vitamin B5 (pantothenic acid)
5. Folic Acid
6. Vitamin B12 (cyanocobalamin)
7. Vitamin H (biotin)
8. Vitamin B2 (riboflavin)

Pharmaceutical


Sulfonamides



Column: Hypersil GOLD aQ, 5 μ m, 150 \times 4.6 mm
 Mobile Phase: A: 0.1% formic acid
 B: Acetonitrile + 0.1% formic acid
 Gradient: 1 0–100% B in 15 min
 Flow Rate: 1.0 mL/min
 Detection: UV @ 270 nm
 Temperature: 30 $^{\circ}$ C

1. Sulfaguanidine
2. Sulfanilamide
3. Sulfathiazole
4. Sulfamerazine
5. Sulfamonomethoxine
6. Sulfaquinoxaline

Xanthines




Column: Hypersil GOLD aQ, 5 μ m, 150 \times 4.6 mm
 Mobile Phase: A: 50 mM monosodium phosphate pH 2.5
 B: Methanol
 Gradient: 1–100% B in 10 min
 Flow Rate: 1 mL/min
 Detection: UV @ 254 nm
 Temperature: 30 $^{\circ}$ C

1. Hypoxanthine
2. Xanthine
3. Theobromine
4. Theophylline
5. Caffeine

Biochemical

PTH amino acids




Column: Hypersil GOLD, 5 μ m, 150 \times 4.6 mm
 Mobile Phase: A: 0.1% tetrahydrofuran + 0.015% triethylamine in water
 B: 0.1% tetrahydrofuran + 0.015% triethylamine in acetonitrile
 Gradient: Time (min) % B
 0 17
 2 20
 7 35
 20 35
 Flow Rate: 1 mL/min
 Detection: UV @ 269 nm
 Temperature: 25 $^{\circ}$ C

1. Serine
2. Asparagine
3. Aspartic acid
4. Glutamic acid
5. Alanine
6. Tyrosine
7. Methionine
8. Tryptophan
9. Phenylalanine
10. Leucine

Hypersil GOLD PFP

Unique selectivity with perfluorinated columns



Particle Size 1.9 μm , 3 μm , 5 μm **USP** L43

Alternative Selectivity to C18 with Excellent Peak Shape and Sensitivity

Hypersil GOLD PFP (pentafluorophenyl) columns build on the performance of Hypersil GOLD silica by providing excellent peak shapes while also offering alternative selectivity in reversed phase chromatography compared to alkyl chain phases. The Hypersil GOLD PFP manufacturing process provides improvements in speed of analysis, peak shape and sensitivity over other fluorinated phases.

Extra Retention for Halogenated Species

Introduction of fluorine groups into the stationary phase causes significant changes in solute-stationary phase interactions. This can lead to extra retention and selectivity for positional isomers of halogenated compounds.

- Analyzing difficult to resolve mixtures of halogenated compounds
- Non-halogenated polar aromatic compounds
- Analysis of complex taxane samples

Unique Selectivity for Non-Halogenated Polar Compounds


Hypersil GOLD PFP Columns are also well suited to the selective analysis of non-halogenated compounds, in particular polar compounds containing hydroxyl, carboxyl, nitro, or other polar groups. High selectivity is often most apparent when the functional groups are located on an aromatic or other rigid ring system.



Hypersil GOLD PFP columns are particularly suited to the analysis of compounds containing substituted aromatic rings. This is because the fluorine atoms around the phenyl ring enhance pi-pi interactions increasing retention and selectivity.

Pharmaceutical


Taxanes



Column: Hypersil GOLD PFP, 5 μ m, 150 \times 4.6 mm
 Mobile Phase: A: Water
 B: Methanol/acetonitrile (7:93)
 Gradient: Time (min) % B
 0 35
 7 35
 25 58
 Flow Rate: 1.5 mL/min
 Detection: UV @ 220 nm

1. 10-Deacetyl baccatin
2. Baccatin III
3. 10-Deacetyl-7-xylosyl taxol B
4. Taxinine M
5. 10-Deacetyl-7-xylosyl taxol
6. 10-Deacetyl taxol
7. 10-Deacetyl-7-xylosyl taxol C
8. 7-Xylosyl taxol
9. Cephalomanine
10. 10-Deacetyl-7epitaxol
11. Paclitaxol
12. Taxol C
13. 7-Epitaxol

Fluorinated nucleic bases




Column: Hypersil GOLD PFP, 5 μ m, 150 \times 4.6 mm
 Mobile Phase: Water + 0.1% tetrahydrofuran
 Flow Rate: 1.0 mL/min
 Temperature: 30 $^{\circ}$ C
 Detection: UV @ 220 nm

1. Fluorocytosine
2. Fluorouracil

Environmental


Banned aromatic amines



Column: Hypersil GOLD PFP, 1.9 μ m, 50 \times 2.1 mm
 Mobile Phase: A: 25 mM ammonium acetate pH 5.0
 B: Acetonitrile
 Gradient: 10–100% B in 3 mins
 Flow Rate: 0.5 mL/min
 Temperature: 40 $^{\circ}$ C
 Detection: UV @ 254 nm (2 μ L flow cell)
 Injection Volume: 0.5 μ L

1. 2,4-Diaminotoluene
2. o-Toluidine
3. 4,4-Oxydianiline
4. 2-Methoxy-5-Methylaniline
5. 2,4,5-Trimethylaniline
6. 4,4-Methylene-bis(2-chloroaniline)
7. Impurity from Analyte No. 6

Phenolic positional isomers






Column: Hypersil GOLD PFP, 5 μ m, 150 \times 4.6 mm
 Mobile Phase: A: Water + 0.1% formic acid
 B: Acetonitrile + 1.0% formic acid
 Gradient: 15–45% B in 20 mins
 Flow Rate: 1.5 mL/min
 Temperature: 25 $^{\circ}$ C
 Detection: UV @ 270 nm
 Injection Volume: 5 μ m

1. 3,4-Dimethoxyphenol
2. 2,6-Dimethoxyphenol
3. 2,6-Difluorophenol
4. 3,5-Dimethoxyphenol
5. 2,4-Difluorophenol
6. 2,3-Difluorophenol
7. 3,4-Difluorophenol
8. 3,5-Dimethoxyphenol
9. 2,6-Dimethoxyphenol
10. 2,6-Dichlorophenol
11. 4-Chloro-3-Methylphenol
12. 3,4-Dichlorophenol
13. 4-Chloro-2-Methylphenol
14. 3,5-Dichlorophenol

Hypersil GOLD CN

Cyano columns for reversed and normal phase separations



Particle Size 1.9 μm , 3 μm , 5 μm **USP** L10


- Steroids and polyphenols in reversed phase
- Surfactants and other polar species in normal phase

Alternative Selectivity with Lower Hydrophobicity than C18

Hypersil GOLD CN columns offer alternative selectivity in reversed phase chromatography with lower hydrophobicity compared to C18 alkyl chain phases. Hypersil GOLD CN columns can also be used in normal phase chromatography, where they offer less retention and different selectivity compared to silica columns.

Pharmaceutical

Penicillins




Column: Hypersil GOLD CN, 5 μm , 150 \times 4.6 mm
Mobile Phase: A: 10 mM potassium phosphate pH3
B: Acetonitrile
Gradient:

Time (mins)	% B
0	0
1	10
8	70

Flow Rate: 1.25 mL/min
Temperature: 25 $^{\circ}\text{C}$
Detection: UV @ 220 nm

1. N-acetyl Penicillamine
2. Ampicillin
3,4,5. Impurities from Penicillin G
6. Penicillin G

TB Drugs



Column: Hypersil GOLD CN, 5 μm , 150 \times 4.6 mm
Mobile Phase: A: 20 mM ammonium formate pH3
B: Acetonitrile
Gradient:


Time (mins)	% B
0	0
15	20

Flow Rate: 1.0 mL/min
Temperature: 25 $^{\circ}\text{C}$
Detection: UV @ 254 nm

1. Isoniazid
2. Pyrazinamide


Toxicology

Steroids



Column:	Hypersil GOLD CN, 5 μ m, 150 \times 4.6 mm	1. Hydrocortisone
Mobile Phase:	A: Water	2. Cortisone
	B: Acetonitrile	3. Corticosterone
Gradient:	Time (mins)	% B
	0	10
	15	50
Flow Rate:	1.5 mL/min	4. 11- α Hydroxprogesterone
Temperature:	25 $^{\circ}$ C	5. 17- α Hydroxprogesterone
Detection:	UV @ 254 nm	6. Progesterone
		7. Deoxycorticosterone


Organic acids



Column:	Hypersil GOLD CN, 5 μ m, 150 \times 4.6 mm	1. 4-Fluorobenzoic
Mobile Phase:	A: 25 mM potassium phosphate pH2	2. o-Toluic Acid
	B: Methanol	3. p-Toluic Acid
Isocratic:	95% A: 5% B	4. 2,4,6-Trimethylbenzoic Acid
Flow Rate:	1.5 mL/min	5. 2,5-Dimethylbenzoic Acid
Temperature:	25 $^{\circ}$ C	
Detection:	UV @ 230 nm	

Hypersil GOLD Phenyl

Excellent retention and unique selectivity for aromatic analytes



- Analyte mixtures with varying polarity and aromaticity
- Where alternative selectivity to C18 is required

Alternative Selectivity for Aromatic and Moderately Polar Analytes

Hypersil GOLD Phenyl reversed phase HPLC columns exhibit alternative selectivity to alkyl chain columns, particularly for aromatic and moderately polar analytes.

Enhanced Pi-Pi Interactions with Aromatics


Many phenyl phases use a propyl (C3) linker between the silica and the phenyl ring. The Hypersil GOLD Phenyl bonded phase contains a butyl (C4) linker which allows for superior alignment of the phenyl ring with aromatic molecules, enhancing pi-pi interactions and therefore their retention.

Moderate Hydrophobicity

The C4 linker also provides the stationary phase with moderate hydrophobicity, making it ideal for the separation of analyte mixtures with varying polarity and aromaticity.

Pharmaceutical


Antibacterials



Column: Hypersil GOLD Phenyl, 5 μ m,
150 \times 4,6 mm
Mobile Phase: A: 20 mM potassium phosphate pH 2.5
B: Acetonitrile
Gradient: 20–50% B in 15 mins
Flow Rate: 1 mL/min
Temperature: 30 $^{\circ}$ C
Injection Volume: 5 μ L
Detection: UV @ 225 nm

1. Carbadox
2. Thiampenicol
3. Furazolidone
4. Oxolinic Acid
5. Sulfadimethoxine
6. Sulfaquinoxaline
7. Nalidixic Acid
8. Piromidic Acid


Antacids



Column: Hypersil GOLD Phenyl, 5 μ m,
150 \times 4,6 mm
Mobile Phase: 20 mM potassium phosphate pH 7.0/
acetonitrile (80/20)
Flow Rate: 1 mL/min
Temperature: 25 $^{\circ}$ C
Injection Volume: 5 μ L
Detection: UV @ 254 nm

1. Famotidine
2. Cimetidine
3. Ranitidine


Veterinary drug coccidiostats



Column: Hypersil GOLD Phenyl, 5 μ m,
150 \times 4,6 mm
Mobile Phase: A: Water
B: Methanol
Gradient: 40–70% B in 10 mins
Flow Rate: 1 mL/min
Temperature: 25 $^{\circ}$ C
Injection Volume: 5 μ L
Detection: UV @ 260 nm

1. 4-amino-3,5-dinitrobenzamide
2. Zoalene (3,5-nitro-o-toluamide)
3. Nitromid (3,5-dinitrobenzamide)
4. Ethopabate

Antidepressants




Column: Hypersil GOLD Phenyl, 1.9 μ m,
50 \times 2.1 mm
Mobile Phase: A: 0.1% formic acid
B: 0.1% formic acid in acetonitrile
Gradient: 10–60% B in 3.4 min
60–90% B in 0.24 min
Flow Rate: 0.5 mL/min
Temperature: 60 $^{\circ}$ C
Injection Volume: 0.7 μ L
Detection: UV @ 225 nm and 254 nm

1. Uracil
2. Acetaminophen
3. p-Hydroxybenzoic acid
4. o-Hydroxybenzoic acid
5. Oxazepam
6. Diazepam
7. Di-isopropyl phthalate
8. Di-n-propyl phthalate

Hypersil GOLD Amino

Highly versatile aminopropyl stationary phase



- Retains anions and organic acids in weak anion exchange
- Excellent for carbohydrate analysis in HILIC

Excellent Chromatographic Properties in Four Modes: Weak Anion Exchange, Reversed Phase, HILIC and Normal Phase


Hypersil GOLD Amino columns can be used with common buffers and an organic modifier as a weak ion exchange material for the analysis of anions and organic acids. When used under normal phase conditions, Hypersil GOLD Amino columns offer an alternative selectivity to silica. Hypersil GOLD Amino columns excel for carbohydrate analysis when used in HILIC mode.

Outstanding Peak Shape and Sensitivity

Based on the same highly pure silica backbone, Hypersil GOLD Amino columns offer improved peak shape over type A silica columns. For high speed, high efficiency separations, Hypersil GOLD Amino columns are available with 1.9 μm particle size.

Food Safety


Sugars



Column: Hypersil GOLD Amino, 5 μm,
150 × 4.6 mm
Mobile Phase: Acetonitrile/water (80:20)
Flow Rate: 1.2 mL/min
Temperature: 35 °C
Detection: RI
Injection Volume: 20 μL

1. Fructose
2. Glucose
3. Sucrose
4. Maltose
5. Lactose

Sorbitol




Column: Hypersil GOLD Amino, 5 μm,
150 × 4.6 mm
Mobile Phase: Acetonitrile/water (80:20)
Flow Rate: 1.2 mL/min
Temperature: 35 °C
Detection: RI
Injection Volume: 20 μL

1. Sorbitol

Hypersil GOLD AX

Separation of anionic species and polar molecules



- Smaller proteins and peptides
- Anionic species
- Polar molecules


Weak Anion Exchange Phase

Hypersil GOLD AX columns utilise a novel polymeric amine ligand bonded to highly pure base deactivated silica. The silica substrate brings higher efficiency than polymer based ion exchange columns.

Suitable for HILIC

Hypersil GOLD AX columns are particularly suited to the analysis of polar compounds in HILIC applications. For high speed, high efficiency separations, Hypersil GOLD AX columns are available with 1.9 μm particle size.


Biopharma Monophosphates



Column: Hypersil GOLD AX, 5 μm, 150 × 4.6 mm
Mobile Phase: Aqueous phosphate buffer (50 mM, pH 3)
Flow Rate: 1.0 mL/min
Temperature: 40 °C
Detection: UV @ 254 nm
Injection Volume: 10 μL

1. Uracil
2. Cytidine-5'-monophosphate
3. Adenosine-5'-monophosphate
4. Guanosine-5'-monophosphate

Food Safety Vitamin C




Column: Hypersil GOLD AX, 5 μm, 100 × 4.6 mm
Mobile Phase: 100 mM ammonium acetate pH 6.8/
acetonitrile (30:70)
Flow Rate: 0.5 mL/min
Temperature: 30 °C
Detection: UV @ 240 nm
Injection Volume: 50 μL

1. Vitamin C

Hypersil GOLD SAX

Quaternary amine strong anion exchange column



- Smaller organic molecules
- Ionic species

High Stability to Aqueous Mobile Phase


The Hypersil GOLD SAX stationary phase utilises a highly stable quaternary amine strong anion exchange ligand bonded to highly pure silica. Hypersil GOLD SAX columns are suited to the analysis of smaller organic molecules such as nucleotides and organic acids using aqueous and low pH mobile phases.

Outstanding Peak Shape and Sensitivity

Based on the same highly pure silica backbone, Hypersil GOLD SAX columns offer improved peak shape over type A silica columns. For high speed, high efficiency separations, Hypersil GOLD SAX columns are available with 1.9 μm particle size.

Biopharma

Monophosphates




Column: Hypersil GOLD SAX, 5 μm ,
150 \times 4.6 mm
Mobile Phase: 50 mM phosphate buffer
pH 3.0
Flow Rate: 1.0 mL/min
Detection: UV @ 254 nm
Column Temperature: 40 $^{\circ}\text{C}$
Injection Volume: 10 μL

1. Uracil
2. Cytidine-5'-monophosphate
3. Adenosine-5'-monophosphate
4. Guanosine-5'-monophosphate

Hypersil GOLD Silica

Excellent peak shape in normal phase chromatography




- Steroids in normal phase
- Polar analytes in HILIC

Outstanding Peak Shape and Sensitivity

Unbonded, highly pure base deactivated silica media that is the backbone of the Hypersil GOLD range of columns. Hypersil GOLD Silica columns are a powerful and efficient tool for the chromatography of non-polar and moderately polar organic compounds by normal phase chromatography. For high speed, high efficiency separations, Hypersil GOLD Silica columns are available with 1.9 μm particle size.

Forensics

Steroids




Column: Hypersil GOLD Silica, 5 μm ,
150 \times 4.6 mm
Mobile Phase: Hexane/ethanol (19:1)
Flow Rate: 1.5 mL/min
Temperature: 30 $^{\circ}\text{C}$
Detection: UV @ 254 nm
Injection Volume: 5 μL

1. Progesterone
2. 21-Hydroxyprogesterone
-21-acetate
3. 17- α -Hydroxyprogesterone
4. Cortisone
5. 11- α -Hydroxyprogesterone
6. Corticosterone
7. Hydrocortisone

Hypersil GOLD HILIC

Enhanced retention of polar and hydrophilic analytes



- Polar and hydrophilic compounds
- Carbohydrates
- Enhanced sensitivity in MS

Enhanced Retention of Polar and Hydrophilic Analytes


Hydrophilic interaction liquid chromatography (HILIC) is an increasingly popular technique offering complementary selectivity to reversed-phase. With the ability to retain highly polar and hydrophilic compounds, Hypersil GOLD HILIC columns have been developed to aid the analysis of compounds that are traditionally difficult to retain using conventional C18 columns. In HILIC, by incorporating water in the highly organic mobile phase, an adsorbed water-rich layer is formed on the polar stationary phase surface into which analyte molecules partition. Retention is governed by dipole-dipole interactions and hydrogen bonding mechanisms.

Improved Sensitivity with MS Detection

The highly organic mobile phases containing low salt levels used for HILIC, make Hypersil GOLD HILIC columns ideal for use with electrospray mass spectroscopy.

Food Safety

Water soluble vitamins




Column: Hypersil GOLD HILIC, 5 μm ,
150 \times 4.6 mm
Mobile Phase: Water/acetonitrile (10:90)
+ 0.1% formic acid
Flow Rate: 1.0 mL/min
Temperature: ambient
Detection: UV @ 205, 230 & 260 nm
Injection Volume: 10 μL

1. Thiamine
2. Nicotinic acid
3. Nicotinamide
4. Pyridoxine
5. Riboflavin
6. PABA.

Chemical

Urea



Column: Hypersil GOLD HILIC, 5 μm ,
150 \times 4.6 mm
Mobile Phase: Water/acetonitrile (10:90)
+ 0.1% formic acid
Flow Rate: 0.6 mL/min
Temperature: 30 $^{\circ}\text{C}$
Detection: +ESI
Injection Volume: 1 μL (made up in mobile phase)

1. Urea

Hypersil GOLD 1.9 μm


Small particles to improve speed and efficiency

The Power of 1.9 μm Particles

1.9 μm particles give higher efficiency than 3 μm or 5 μm particles and this efficiency is delivered over a greater range of optimum linear velocity. This makes it possible to operate at higher flow rates without losing performance. Because shorter columns packed with 1.9 μm particles give equivalent efficiency to longer columns packed with 5 μm particles, faster analysis and solvent savings for the chromatographer become a reality.

Three Tips for Method Transfer

1. To maintain an equivalent separation when transferring a method it is important to keep the reduced linear velocity constant between the original and new method.
2. Sub-2 μm based methods are most often transferred to smaller volume columns, so the same injection volume will take up a larger proportion of the new column, possibly leading to band broadening. It is therefore important to scale down the injection volume to match the change in column volume.
3. Geometrical transfer of the gradient requires calculation of the number of column volumes of mobile phase in each segment (time interval) of the gradient in the original method to ensure that the new calculated gradient takes place over the same number of column volumes, for the new column.




HPLC Method Development Calculator
www.thermoscientific.com/crc

Pressure Rating of Hypersil GOLD 1.9 μm Columns

Column Hardware	Pressure Rating
Analytical columns	1250 bar/18,000 psi
Capillary/nano columns	400 bar/6,000 psi
Javelin HTS columns	400 bar/6,000 psi


Transferring a method using these tips can give results as shown below for the separation of Ibuprofen and impurities.



Which 1.9 μm Column?

We offer an extensive range of columns packed with 1.9 μm particles to suit the full variety of application needs. The choice of column will depend upon the requirement of the analysis.

- Speed: choose from 10 mm Javelin HTS, 20, 30 or 50 mm long analytical columns
- Efficiency: choose a longer column (for example 150 or 200 mm)
- Low backpressure: Hypersil GOLD 1.9 μm media is packed into a high pressure column 50 mm long and 4.6 mm internal diameter. Traditionally, a 1.9 μm column is used on UHPLC instruments. However, by producing less backpressure, this new wider column is suitable for users of conventional systems where pressure limits are often in the 6000 psi/400 bar region, ensuring fast chromatography without the need for extensive instrument optimization.




Hypersil GOLD 1.9 μm Javelin HTS Columns for Speed

Hypersil GOLD 1.9 μm Javelin HTS columns take fast LC to the extreme. These short 10 mm columns enable analysis times as fast as 8 seconds to be achieved. The use of ultra-low dead volume, direct connect Javelin hardware also minimizes dispersion.

Toxicology

Nandrolone




Column:	Hypersil GOLD 1.9 μm , 10 \times 2.1 mm	1. Nandrolone (19-Nortestosterone)
Mobile Phase:	Water/acetonitrile 40/60 + 0.1% tetrahydrofuran isocratic	
Flow Rate:	0.4 mL/min	
Temperature:	5 $^{\circ}\text{C}$	
Detection:	254 nm	
Injection Volume:	0.5 μL	

200 mm Column for Efficiency

The 1.9 μm particles used in Hypersil GOLD columns give less backpressure than 1.8 or 1.7 μm , permitting the use of longer columns for greater efficiency.

Environmental

Phenolic pollutants



Column:	Hypersil GOLD 1.9 μm , 200 \times 2.1 mm	1. Phenol 2. 4-nitrophenol 3. 2-nitrophenol 4. 4-chlorophenol 5. 2-chlorophenol 6. 2,6-dimethylphenol 7. 2,4-dimethylphenol 8. 2,4-dichlorophenol 9. Pentachlorophenol
Mobile Phase:	A: Water + 0.1% formic acid B: Acetonitrile + 0.1% formic acid	
Gradient:	5 to 95% B in 24 min	
Flow Rate:	0.6 ml/min	
Temperature:	60 $^{\circ}\text{C}$	
Injection Volume:	1 μL	
Detection:	UV at 270 nm	
Pressure:	606 bar	

System Considerations

With 1.9 μm particles, analyses can be performed with a high linear velocity through the column without loss in performance, provided the LC system is optimized to operate under these conditions. In order to produce fast, efficient chromatography, all system components for the assay should also be considered. Modern ultra high pressure liquid chromatography (UHPLC) instruments, including the Thermo Scientific™ Vanquish™ UHPLC system, will take account of these factors.

There are three major system considerations to remember when using short columns packed with 1.9 μm particles.

1. The system volume (connecting tubing ID and length, injection volume, UV detector flow cell volume) must be minimized
2. The detector time constant and sampling rate need to be carefully selected
3. When running fast gradients pump delay volume needs to be minimal.

Hardware Solutions

Hardware solutions for high throughput screening, capillary and preparative chromatography

Hypersil GOLD columns are available in particle sizes and column designs to meet all separation needs, including improved resolution, enhanced sensitivity, and faster analyses. With particle sizes from 1.9 μm to 12 μm , Hypersil GOLD columns offer chromatographic solutions with consistent separations and performance. Specialized hardware includes:

- Preparative columns
- Thermo Scientific™ Javelin™ HTS direct-connection columns
- Guard columns for column protection

Preparative Columns

Analytical methods may require scale up to preparative sizes to isolate and purify compounds from mixtures. In choosing the best column and packing material for your preparative application, consider:

- Selectivity
- Loadability of the media
- Column dimensions

We have established a strong reputation for the manufacture and supply of high quality preparative columns, designed to give the same levels of performance and reproducibility as our popular analytical columns. Scale up is easiest when starting from an analytical column packed with smaller particle size media offering the same selectivity as the larger particle size preparative media. Hypersil GOLD phases are offered in various sizes to complement lab scale operations and facilitate the scale up to preparative chromatography. Contact us for ordering details on Hypersil GOLD preparative columns.



Columns for High Throughput Screening

Javelin HTS columns are specifically designed for high throughput applications. Using finger tight fittings and low dead volume hardware to minimise band broadening, these columns are ideal for ballistic gradients, providing enough retention and sensitivity for very fast assays. Javelin HTS columns are available in multipacks to provide a cost effective solution.



Javelin HTS Column

Description	Particle Size	Length (mm)	ID (mm)	Part Number
Hypersil GOLD Javelin HTS Column (3/pk)	1.9	10	2.1	11808901
	5	20	4.0	10157354

Column Protection

Extend column lifetime and improve performance

Guard Columns

Drop-in guard cartridges and holders offer convenience, economy, and effective protection for extending analytical column lifetimes. The 10 mm design offers maximum protection with minimal increase in retention. Hypersil GOLD drop-in guard cartridges are provided in packs of 4 each.



UHPLC Filter

Replaceable 0.2 µm Thermo Scientific UHPLC filter cartridges can be used to protect Hypersil GOLD 1.9 µm columns against particulate contamination, extending column lifetime. Its low dead volume design maintains chromatographic performance without degrading peak shape and causes minimal efficiency loss through dispersion. The UHPLC filter adds minimal increase in backpressure and so can be fitted to any length column.



UNIGUARD Guard Cartridge Holder

Description	Length (mm)	ID (mm)	Part Number
UNIGUARD Guard Cartridge Holder	10	1.0	10762567
		2.1	10776714
		3.0	10776714
		4.0/4.6	10602864

Description	ID (mm)	Part Number
UHPLC Filter Holder		10775706
UHPLC Filter Cartridge, 0.2 µm (5/pk)	2.1	10127594
	1.0	11550944

Ordering Information

Hypersil GOLD HPLC Columns

Particle Size (µm)	Description	Length (mm)	ID (mm)	Hypersil GOLD	Hypersil GOLD C8	Hypersil GOLD C4			
1.9	UHPLC Column	20	2.1	10464613	10658895	–			
			3.0	10305063	–	–			
		30	2.1	10049142	10044264	–			
			1.0	10020213	–	–			
			2.1	10474783	10664945	10127614			
			3.0	10345203	10556845	–			
			4.6	10169344	10096024	–			
			100	1.0	10715004	10329884	–		
		150	2.1	10734604	10074694	10691396			
			3.0	10582465	10442195	–			
			200	2.1	10630204	10234864	10391585		
		3	Drop-in Guard (4/pk)	10	1.0	10608705	10714817	10147884	
2.1	11518270				11598270	10447534			
3.0	10481365				10075814	10157884			
4.0/4.6	11528270				10487694	10610626			
HPLC Column	30		2.1	10394463	–	–			
			3.0	–	–	–			
			4.6	10764984	–	–			
	50		2.1	10424413	10553035	10159534			
			3.0	10293183	10558585	–			
			4.0	10658825	10715126	–			
			4.6	10089282	10543995	–			
	100		1.0	10281225	10027464	–			
			2.1	10632864	10018912	10037254			
			3.0	10304923	10619395	10027504			
			4.0	10794546	10202105	–			
			4.6	10384323	10774414	10732826			
			150	1.0	10704556	–	10199474		
	150		2.1	10121223	10078952	10599735			
			3.0	10131383	10312745	10312355			
			4.0	10456844	10734766	–			
			4.6	10610564	10736144	10261325			
			5	Drop-in Guard (4/pk)	10	2.1	11538270	10535445	10076694
						3.0	10544085	10259644	10514475
	4.0/4.6					11548270	11508280	10670046	
	HPLC Column			30	2.1	10475163	–	–	
					3.0	10668825	–	–	
					4.6	10008862	–	–	
				50	2.1	10611334	10293413	10609585	
					3.0	10191523	10764356	–	
					4.6	10611524	10690564	10699385	
				100	2.1	10140743	10503235	10382865	
					3.0	10384723	10261915	10322355	
					4.6	10121033	10690754	10699755	
			150	2.1	10131033	10203663	10609765		
				3.0	10795316	10108874	–		
				4.0	10313065	10149494	–		
4.6		10501695		10755783	10486854				
250		2.1		10009102	10765783	10067684			
		3.0		10352775	10619405	10548415			
		4.0	10506365	10343465	–				
		4.6	10543035	10243453	10312635				

Hypersil GOLD HPLC Columns

Particle Size (µm)	Description	Length (mm)	ID (mm)	Hypersil GOLD aQ	Hypersil GOLD PFP	Hypersil GOLD CN		
1.9	UHPLC Column	20	2.1	10375203	10670774	–		
			3.0	–	–	–		
		30	2.1	10324223	10243553	–		
			1.0	10101723	10465753	–		
			2.1	10572635	10661134	10319984		
			3.0	10111723	10435943	–		
	50	1.0	10621344	10523435	–			
		2.1	10582635	10583025	10565885			
	3	Drop-in Guard (4/pk)	10	1.0	10360905	10486764	10269954	
				2.1	11548280	11528290	11588290	
				3.0	11558280	10734657	10035844	
				4.0/4.6	11568280	10270055	11598290	
HPLC Column		30	2.1	10715793	10485533	10343993		
			3.0	10466664	10260995	–		
5	Drop-in Guard (4/pk)	10	2.1	10783497	10619285	10361105		
			3.0	10370905	10535825	10341215		
			4.0/4.6	11578280	11538290	11508300		
			HPLC Column	30	2.1	–	–	–
					3.0	–	–	–
					4.6	–	–	–
	50	2.1	10776153	10131223	10600954			
		3.0	10508405	10640536	–			
		4.6	10662474	–	10652484			
	100	2.1	10140593	10620564	10662484			
		3.0	10557435	10704756	10231965			
		4.6	10150593	10621524	–			
	150	2.1	2.1	10718063	10724414	10746144		
			3.0	10497414	10220815	–		
			4.0	–	10418574	–		
		4.6	4.6	10019382	10068902	10213663		
			2.1	10161023	10078902	10474983		
			3.0	–	–	10027154		
250	4.0	10056134	10159194	10538585				
	4.6	10582065	10141033	10273703				

Hypersil GOLD HPLC Columns

Particle Size (µm)	Description	Length (mm)	ID (mm)	Hypersil GOLD Phenyl	Hypersil GOLD Amino	Hypersil GOLD AX	
1.9	UHPLC Column	20	2.1	–	–	–	
			30	1.0	–	–	–
			2.1	–	–	–	
		50	1.0	–	–	–	
			2.1	10422015	10793097	10449644	
			3.0	–	–	–	
			4.6	–	–	–	
			100	1.0	–	–	–
				2.1	10207834	10177994	10601406
		150	2.1	10431825	10670246	10229564	
			200	2.1	10595125	10732347	10611406
		3	Drop-in Guard (4/pk)	10	1.0	10459054	10036304
2.1	10380915				10698515	10016414	
3.0	10279654				10478664	10742137	
4.0/4.6	10026124				10056114	10167734	
HPLC Column	30		2.1	–	10681196	10077194	
			3.0	–	–	–	
			4.6	–	–	–	
	50		2.1	10208264	10187654	10037824	
			3.0	–	–	–	
			4.0	–	–	–	
			4.6	–	–	–	
	100		1.0	–	–	–	
			2.1	10056364	10742527	10097834	
			3.0	10437324	10689755	–	
			4.0	–	–	–	
			4.6	10690256	10289364	10467254	
150			1.0	10794627	10211095	10087344	
			2.1	–	10600826	10763606	
			3.0	10446984	10169564	10139184	
	4.0		–	–	–		
5	Drop-in Guard (4/pk)		10	2.1	10025644	10659465	10469644
				3.0	10525435	10525615	10536185
				4.0/4.6	10351445	10259174	10734817
	HPLC Column		30	2.1	–	–	–
		3.0		–	–	–	
		4.6		–	–	–	
		50	2.1	10675335	10427324	–	
			3.0	–	–	–	
			4.6	10762907	10157554	10547245	
		100	2.1	10753477	10391295	10302455	
			3.0	10127604	10169524	–	
			4.6	–	–	10231385	
		150	2.1	–	–	10567455	
			3.0	–	–	–	
			4.0	–	–	–	
			4.6	10014574	10006074	10077364	
			250	2.1	–	10554105	10332885
				3.0	10594095	10517455	10734366
		4.0		10005784	10322975	10648995	
		4.6		10442595	10752717	–	

Hypersil GOLD HPLC Columns

Particle Size (µm)	Description	Length (mm)	ID (mm)	Hypersil GOLD SAX	Hypersil GOLD Silica	Hypersil GOLD HILIC	
1.9	UHPLC Column	20	2.1	–	–	–	
			30	1.0	–	–	–
			2.1	–	–	–	
		50	1.0	–	–	–	
			2.1	–	10722157	11527931	
			3.0	–	–	–	
			4.6	–	–	–	
			100	1.0	–	–	–
		100	2.1	10259514	10065884	11567931	
			3.0	–	–	–	
			150	2.1	10754257	10341165	10029259
200	2.1	–	10498674	–			
3	Drop-in Guard (4/pk)	10	1.0	–	10468664	11738056	
			2.1	10620056	10371105	11748056	
			3.0	10056224	10629865	11758056	
			4.0/4.6	10609685	10289124	11768056	
			HPLC Column	30	2.1	–	10524125
	3.0	–	–		–		
	4.6	–	10076744		–		
	50	2.1	10362645	10187524	11597941		
		3.0	–	–	–		
		4.0	–	–	–		
		4.6	–	–	–		
		100	1.0	–	–	–	
		2.1	10609025	10066264	11547951		
		3.0	10487414	10177954	11557951		
	4.0	–	–	–			
	150	4.6	10569155	–	11577951		
		1.0	–	10763887	11587951		
		2.1	–	10764447	11597951		
		3.0	10540126	10127624	11507961		
		4.0	–	–	–		
		4.6	10392475	10157894	11527961		
		5	Drop-in Guard (4/pk)	10	2.1	10702337	10555425
	3.0	10703697			10702147	11718056	
	4.0/4.6	10006844			10648905	11728056	
	5	HPLC Column	30	2.1	–	–	–
				3.0	–	–	–
				4.6	–	–	–
			50	2.1	10417454	10311535	11537971
				3.0	–	–	–
4.6				–	–	11567971	
100			2.1	10352525	10311125	11587971	
			3.0	10251185	10418284	11597971	
			4.6	10477414	10290385	–	
150		2.1	10598195	10167474	–		
		3.0	–	–	–		
		4.0	–	–	–		
		4.6	10629755	10437534	11567981		
		250	2.1	10149084	10556955	11587981	
		3.0	10417624	–	11597981		
250		4.0	10333065	10774077	11507991		
		4.6	10703786	10360835	11517991		


Resources

for Chromatographers

Thermo Scientific Chromatography Columns and Consumables Catalog

This extensive catalog offers 450 pages of proven chromatography tools and product selection guides. Available online, with a robust search tool and optimized for your iPad®.

Visit www.thermofisher.com/catalog



Chromatography Resource Center

Our web-based resource center provides technical support, applications, technical tips and literature to help move your separations forward.

Visit www.thermofisher.com/sbe



Chromexpert

A dedicated team of technical support specialists and application chemists available to help you with product selection and assistance when using your chromatography consumables. Our experts have access to the latest chromatography technology and will act as your trusted advisors. For application and technical support.

Visit www.thermofisher.com/chromexpert



AppsLab Library

Thermo Scientific™ AppsLab Library of Analytical Applications provides more than 1300 detailed application examples for the columns listed in the 2016–2017 Chromatography Columns and Consumables Catalog. Search, filter and download complete methods to optimize your separation or implement validated methods using Thermo Scientific™ Dionex™ Chromeleon™ Chromatography Data System (CDS) software. AppsLab Library makes our global application expertise accessible to you—online and downloadable.

Visit www.thermofisher.com/AppsLab



thermo scientific



For more information, contact Fisher Scientific



Austria: +43 (0)800 20 88 40 **Belgium:** +32 (0)056 260 260 **Denmark:** +45 70 27 99 20
Germany: +49 (0)2304 932-5 **Ireland:** +353 (0)1 885 5854 **Italy:** +39 02 950 59 478
Finland: +358 (0)9 8027 6280 **France:** +33 (0)3 88 67 14 14 **Netherlands:** +31 (0)20 487 70 00
Norway: +47 22 95 59 59 **Portugal:** +351 21 425 33 50 **Spain:** +34 902 239 303
Sweden: +46 31 352 32 00 **Switzerland:** +41 (0)56 618 41 11 **UK:** +44 (0)1509 555 500

ThermoFisher
SCIENTIFIC

For Research Use Only. Not for use in diagnostic procedures. © 2017 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. **TG21486-EN 0217S**